Qlip control samples programme

Qlip calibration and control samples have been specifically developed for the calibration and control of infrared instruments.

Registration form 2022



| Company | Contact person (1st) | |
|------------------|-----------------------------------|--|
| Delivery | Contact person (2 nd) | |
| address | E-mail (1st) | |
| (Street address) | E-mail (2 nd) | |
| address) | Phone (for courier) | |
| | Phone (2 nd) | |

- Please indicate the number of samples you would like to receive in the yellow fields and if you want to recieve the samples weekly, 2 weekly (even weeks) or 4 weekly (even weeks).

 Please note that you will be charged according to the rate indicated in the 7th column. The costs are for each individual delivery.

 Each sample contains 75-80 m per bottle only somatic cell count is 35-40 ml per bottle.

 Actual rates available on request

| | | | | | | The yellow fields are available | | |
|-----------------|-----------------------|---|-------------------------------------|---------------------------------------|--------------|---------------------------------|----------|----------|
| Type of product | Description | Certified parameters | Approximate value | Amount of samples in the subscription | Article code | Weekly | 2 Weekly | 4 Weekly |
| Milk | raw milk | fat, protein, lactose, total solids, urea, casein, non-protein-nitrogen, freezing point | 4,4% m/m fat, 3,6% m/m protein | 1-15 samples | SC5300e | | | |
| Milk | cheese milk A | fat, protein, lactose, total solids | 3,7% m/m fat, 3,6% m/m protein | 1-15 samples | SC5140e | | | |
| Milk | cheese milk B | fat, protein, lactose, total solids | 2,7% m/m fat, 3,6% m/m protein | 1-15 samples | SC5180e | | | |
| Milk | cheese milk C | fat, protein, lactose, total solids | 1,5% m/m fat, 3,6% m/m protein | 1-15 samples | SC5220e | | | |
| Milk | semi-skimmed milk | fat, protein | 1,5% m/m fat, 3,4% m/m protein | 1-15 samples | SC5100e | | | |
| Milk | whole milk | fat, protein | 3,5% m/m fat, 3,4% m/m protein | | SC5110e | | | |
| Milk | skimmed milk | fat, protein, lactose, total solids | 0,04% m/m fat, 3,6% m/m protein | | SC5260e | | | |
| Sheep milk | sheep milk | fat, protein, lactose, total solids | 6% m/m fat, 5% m/m protein | | SC5136e | | | |
| Goat milk | goat milk | fat, protein | 4% m/m fat, 3,5% m/m protein | | SC5120e | | | |
| Goat milk | goat milk | fat, protein, total solids | 4% m/m fat, 3,5% m/m protein | | SC5125e | | | |
| Goat milk | goat milk | fat, protein, lactose, total solids | 4% m/m fat, 3,5% m/m protein | 1-15 samples | SC5131e | | | |
| Goat milk | skimmed goat milk | fat, protein, total solids | 0,1% m/m fat, 3,6% m/m protein | | SC5139e | | | |
| Goat cream | goat cream | fat, protein, total solids | 42% m/m fat, 2,0% m/m protein | 1-15 samples | SC5138e | | | |
| Goat whey | goat whey | fat, protein | 0,65% m/m fat, 0,94% m/m protein | 1-15 samples | SC5130e | | | |
| Goat whey | goat whey | fat, protein, total solids | 0,65% m/m fat, 0,94% m/m protein | 1-15 samples | SC5135e | | | |
| Goat whey | goat whey | fat, protein, lactose, total solids | 0,65% m/m fat, 0,94% m/m protein | 1-15 samples | SC5132e | | | |
| Goat whey | centrifuged goat whey | fat, protein | 0,1% m/m fat, 0,7% m/m protein | 1-15 samples | SC5423e | | | |
| Whey | 1st whey | fat, protein, lactose, total solids | 0,24% m/m fat, 0,85% m/m protein | 1-15 samples | SC5020e | | | |
| Whey | 2nd whey | fat, protein, lactose, total solids | 0,15% m/m fat, 0,60% m/m protein | 1-15 samples | SC5060e | | | |
| Whey | centrifuged whey | fat, total solids | 0,03% m/m fat, 5% m/m total solids | 1-15 samples | SC5419e | | | |
| Whey | centrifuged whey | fat, protein, total solids | 0,03% m/m fat, 0,6% m/m protein | | SC5417e | | | |
| Whey | whey concentrate | fat, total solids | 0,12% m/m fat, 17% m/m total solids | 1-15 samples | SC5080e | | | |
| Retentate | retentate | fat, protein, lactose, total solids, non-protein-nitrogen, casein | 0,15% m/m fat, 10,5% m/m protein | | SC5330e | | | |
| Cream | cream A | fat, protein, lactose, total solids, freezing point | 42% m/m fat, 2% m/m protein | | SC5340e | | | |
| Cream | cream B | fat | 35% m/m fat | 1-15 samples | SC5380e | | | |
| Cream | cream C | fat, protein | 21% m/m fat, 2,8% m/m protein | 1-15 samples | SC5390e | | | |
| Cream | cream D | fat | 31% m/m fat | 1-15 samples | SC5405e | | | |
| Cream | cream 10% | fat, protein, lactose, total solids | 10% m/m fat | | SC5468e | | | |
| Milk | milk | somatic cell count low | 200 x 1000 cells/ml | 1-15 samples | SC5191e | | | |
| Milk | milk | somatic cell count high | 400 x 1000 cells/ml | 1-15 samples | SC5190e | | | |

| | | | | | | The yellow fields are available | | |
|-----------------|--------------------|---|-------------------------------------|---------------------------------------|--------------|---------------------------------|----------|----------|
| Type of product | Description | Certified parameters | Approximate value | Amount of samples in the subscription | Article code | Weekly | 2 Weekly | 4 Weekly |
| Milk | raw milk | fat, protein, lactose, total solids, urea, casein, non-protein-nitrogen, freezing point | 4.4% m/m fat. 3.6% m/m protein | | SC5305e | | - '' | - 1 |
| Milk | cheese milk A | fat, protein, lactose, total solids, drea, casein, non-protein-nitrogen, neezing point | 3,7% m/m fat, 3,6% m/m protein | 16-50 samples | SC5145e | | | |
| Milk | cheese milk B | fat, protein, lactose, total solids | 2.7% m/m fat, 3,6% m/m protein | | SC5185e | | | |
| Milk | cheese milk C | fat, protein, lactose, total solids | 1.5% m/m fat, 3.6% m/m protein | 16-50 samples | SC5225e | | | |
| Milk | skimmed milk | fat, protein, lactose, total solids | 0.04% m/m fat. 3.6% m/m protein | 16-50 samples | SC5265e | | | |
| Goat milk | goat milk | fat, protein | 4% m/m fat, 3,5% m/m protein | 16-50 samples | SC5430e | | | |
| Goat milk | goat milk | fat, protein, total solids | 4% m/m fat, 3,5% m/m protein | 16-50 samples | SC5124e | | | |
| Goat whey | goat whey | fat, protein | 0,65% m/m fat, 0,94% m/m protein | 16-50 samples | SC5431e | | | |
| Whey | 1st whey | fat, protein, lactose, total solids | 0,24% m/m fat, 0,85% m/m protein | 16-50 samples | SC5025e | | | |
| Whey | 2nd whey | fat, protein, lactose, total solids | 0,15% m/m fat, 0,60% m/m protein | 16-50 samples | SC5065e | | | |
| Whey | centrifuged whey | fat, total solids | 0,03% m/m fat, 5% m/m total solids | 16-50 samples | SC5420e | | | |
| Whey | whey concentrate | fat, total solids | 0,12% m/m fat, 17% m/m total solids | 16-50 samples | SC5085e | | | |
| Cream | cream A | fat, protein, lactose, total solids, freezing point | 42% m/m fat, 2% m/m protein | 16-50 samples | SC5345e | | | |
| Cream | cream B | fat | 35% m/m fat | 16-50 samples | SC5386e | | | |
| Cream | cream C | fat, protein | 21% m/m fat, 2,8% m/m protein | 16-50 samples | SC5395e | | | |
| Cream | cream D | fat | 31% m/m fat | | SC5385e | | | |
| Cream | cream 48% (5 sets) | fat, protein | 48% m/m fat, 1,7% m/m protein | 5x34 samples | SC5160e | | | |
| Cream | cream 10% | fat | 10% m/m fat | 16-50 samples | SC5348e | | | |

| Optional results you can receive | | | Article code | Yes/No |
|----------------------------------|-------------------|---------------------------------|--------------|--------|
| Milk | skimmed milk | freezing point | SC5268e | |
| Milk | semi-skimmed milk | freezing point | SC5106e | |
| Goat milk | goat milk | | SC5133k | |
| Goat milk | skimmed goat milk | | SC5134k | |
| Goat cream | goat cream | | SC5142k | |
| Whey | | protein | SC5093e | |
| Milk | raw milk | acidity fat (BDI), pH | SC5318e | |
| Milk | raw milk | somatic cell count (Fossomatic) | SC5314e | |

| Remarks: | |
|----------|--|
| | |
| | |
| | |